


Lab Connection Instructions

Carestream



DEXIS / Carestream Connect Instructions:

1. Login to the platform as a Dentist (CS Connect or DEXIS IS Connect)
2. In the toolbar, click 
3. In the Partners window, click Invite a partner or Add a partner.
4. In the Email field, enter the email address of the Laboratory

Use exact email: **scans@smiledesign.cloud**

5. In the Message box, either leave the default text or write your own message.
6. Click Send.

The invitation message is sent to the Laboratory and a Partner entry is created in the list of Partners.

7. When the Laboratory logs in to CS Connect or DEXIS IS Connect, the Partners  button displays a notification as follows: